ULTRACOR[®] TWIRL[®]

Breast Tissue Marker

ENGLISH

INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The ULTRACOR[®] TWIRL[®] Breast Tissue Marker consists of a disposable beveled needle applicator containing a Nitinol radiographic ring wireform. The wireform is intended for long-term radiographic marking of the tissue site. The applicator has a beveled 17G x 10cm needle with 1cm depth marks and a locking plunger. The ring is deployed from the beveled needle tip into the tissue site.

INDICATIONS FOR USE

The ULTRACOR[®] TWIRL[®] Breast Tissue Marker is intended for use to attach to soft breast tissue, including axillary lymph nodes, to radiographically mark the location of the biopsy procedure.

CONTRAINDICATIONS

Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant. The implant is made from a nickel-titanium alloy; if there is a known allergy to nickel, use of the ULTRACOR[®] TWIRL[®] Breast Tissue Marker is not advised.

WARNINGS

- **As with any foreign object implanted into the body, potential adverse reactions are possible. It is the responsibility of the physician to evaluate the risk/benefit prior to the use of this device.**
- **The ULTRACOR[®] TWIRL[®] Breast Tissue Marker is not recommended for use in patients with breast implants.**
- **This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.**
- **Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**

PRECAUTIONS

- This product should only be used by a physician who is completely familiar with the indications, contraindications, limitations and typical findings and possible side effects of breast tissue marker placement.
- Do not use if needle is bent and/or tip is damaged.
- Use caution when handling the device to prevent premature deployment of the breast tissue marker.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

POTENTIAL COMPLICATIONS

Potential complications of breast tissue marker placement may include, but are not limited to: hematoma, hemorrhage, infection, allergic reaction, marker migration, misdiagnosis, lymphedema, adjacent tissue injury, pain, and nerve injury.

HOW SUPPLIED

The ULTRACOR[®] TWIRL[®] Breast Tissue Marker is provided sterile, unless the package has been opened or damaged. For single use only.

EQUIPMENT REQUIRED


The following equipment is required to place the breast tissue marker:

- Appropriate imaging modality and accessories
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

DIRECTIONS FOR USE (Refer to Figure 1)


1. Prepare the site as required. Adequate anesthesia should be administered, as required.
2. Inspect the packaging to insure that package integrity has not been compromised. The product is sterile unless the seal is broken.
3. Using standard aseptic technique, remove the ULTRACOR[®] TWIRL[®] Breast Tissue Marker device from the package. **REMOVE THE TIP PROTECTOR (A)** in Figure 1. Inspect the needle tip for signs of damage.
4. Using the 1cm depth markings on the needle (B), insert the introducer needle directing the needle tip to the target site. Note that the BARD[®] logo (D) is aligned with the needle bevel.
5. Confirm needle placement with the appropriate imaging technique. If necessary, reposition the needle and reconfirm placement.
6. Release the plunger safety lock by turning the plunger (E) clockwise 90 degrees from the locked position (F) until the lock symbol is vertically aligned with the logo and the red indicator marks on the applicator and plunger are aligned in the unlocked position (G).
7. Deploy the Marker (H) by advancing the plunger until it contacts the barrel of the applicator (C) while maintaining the applicator stable in position.

NOTE: The marker is designed to deploy proximal to the tip of the needle as illustrated below:

8. Remove the applicator.
9. Dispose of the applicator properly.
10. Confirm final marker position using the appropriate imaging technique.

FIGURE 1
ULTRACOR® TWIRL® Breast Tissue Marker

MRI Safety Information

Non-clinical testing demonstrated that the ULTRACOR® TWIRL® Breast Tissue Marker is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla only
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Mode)

Under the scan conditions defined above, the ULTRACOR® TWIRL® Breast Tissue Marker is expected to produce a maximum temperature rise of 2.3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 10 mm from the ULTRACOR® TWIRL® Breast Tissue Marker when imaged using a gradient echo pulse sequence and a 3 Tesla MRI system.

The delivery system (applicator) for the ULTRACOR® TWIRL® Breast Tissue Marker has not been tested for MR compatibility.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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**Copyright © 2018 by C. R. Bard Inc. All rights reserved.**
Copyright © 2018 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2018 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2018 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2018, C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2018 C. R. Bard, Inc. Alle rechten voorbehouden.